


# Full face protection shield



- 1 Flap\*
- 2 forehead strip\*
- 3 elastic band

\*PET sheet

## Variants:

- A) Full size ("L") or shortened ("S") in low weight (0.3 mm) or standard (0.5 mm) variants
- B) combined with protection mantle
- C) with anti-slip forehead stripe to wear over Tyvek suit



## To sterilize/disinfect:




- remove elastic band (to be incinerated)
- apply antiseptic liquid or UV
- assemble with new elastic band

- ✓ Supplier: Skoltech, Moscow
- ✓ 2000...8000 sets/day\*
- ✓ Only local materials
- ✓ Reusable sterilizable
- ✓ Changeable elastic bands in the kit
- ✓ Low cost
- ✓ Extra light

\* Limited free supply is available by request from hospitals

# How to assemble? (2-3 min)

Apply forehead stripe to flap and insert elastic band into the side slot



Size and mass:

 L size Light — 40 g  
L size Standard — 53g

